A large iceberg floats in the ocean under a sunset sky. The top part of the iceberg is above the water, while a much larger, rectangular part is submerged below the surface. The water is dark blue, and the sky is a mix of orange, yellow, and blue.

# Is contamination of bronchoscopes really a problem?



# Is contamination of bronchoscopes really a problem?

## What we know

**210**

More than 210 patients contaminated by bronchoscopes in US from 2014-2016 (1)

**183**

183 medical device reports to FDA on contaminated bronchoscopes in 2016 (1)

**60**

More than 60 outbreak articles published on contaminated bronchoscopes (2)

### Top 1 Hazard

Contamination and inadequate cleaning of endoscopes on ECRI annual list for "Top 10 Health Technology Hazards" past 8 years (3)


**250**

Between 2012 and the spring of 2015, endoscopes caused at least 250 life-threatening infections worldwide, including infections with the superbug carbapenem-resistant enterobacteriaceae, according to the results of an investigation conducted by a U.S. Senate committee (4)

**More nosocomial infection and pseudo-infection outbreaks have been linked to contaminated endoscopes than to any other medical device. (5)**


# High Level Disinfection does not methodically clean bronchoscopes

Cleaning and High Level Disinfection eliminate somewhere between 6 and 12 logs of microorganisms, but endoscopes can potentially contain 10 log bioburden. Thus even after cleaning and highlevel disinfection, scopes can potentially still have 4 logs left or as many as 10,000 organisms left before essentially next patient use. (6)



## High increase of bronchoscope medical device reports to FDA

No of MDR reports on infection or device contamination associated with reprocessed flexible bronchoscopes (1)



Increased focus on cross-contaminated endoscopes from US authorities over the recent years has led to a high increase in filing of Medical Device Reports (MDR) on bronchoscopes.

The filed reports mention infection or device contamination associated with reprocessed flexible bronchoscopes.

## Flexible endoscopes continuously on ECRI Top 10 Health Technology Hazards list

ECRI top 10 health technology hazards year 2010-2017 (3)


2017	<p><b>2.</b> Inadequate Cleaning of Complex Reusable instruments Can Lead to Infections</p> <p><b>10.</b> Device Failures Caused by Cleaning Products and Practices</p>
2016	<p><b>1. Inadequate Cleaning of Flexible Endoscopes before Disinfection Can Spread Deadly Pathogens</b></p>
2015	<p><b>4.</b> Inadequate Reprocessing of Endoscopes and Surgical Instruments</p>
2014	<p><b>6.</b> Inadequate Reprocessing of Endoscopes and Surgical Instruments</p>
2013	<p><b>8.</b> Inadequate Reprocessing of Endoscopes and Surgical Instruments</p>
2012	<p><b>4.</b> Cross-contamination from Flexible Endoscopes</p>
2011	<p><b>3.</b> Cross-Contamination from Flexible Endoscopes</p>
2010	<p><b>1.</b> Cross-Contamination from Flexible Endoscopes</p>

For the past 8 years cross-contamination of endoscopes has been recognized as a major patient safety issue by the ECRI organization.

ECRI Institute is a nonprofit US organization focusing on improving patient safety.

# Since 2015 US authorities have increased cross-contamination considerably

# their focus on endoscope



# Reprocessing of an endoscope is highly complex with more than 100 steps

Flexible bronchoscopes are difficult to clean and disinfect due to the long and narrow channel. More healthcare-associated outbreaks have been linked to contaminated endoscopes than to any other medical device. (5,12)


Failure in compliance with scientifically-based reprocessing guidelines has led to numerous outbreaks. However, the persistence of contamination on endoscopes, even after adequate reprocessing, is well documented. (8,5,13,14)

In the March 2015 reprocessing guideline, FDA identifies bronchoscopes as being part of a subset of devices that poses a greater likelihood of microbial transmission and represents a high risk of infection if not adequately reprocessed. (15)

With more the than 100 steps for reprocessing each endoscope after use, adherence to new complex guidelines (AAMI/ARON/SGNA) is both costly, time-consuming, and occupying various resources. Cost associated with the reprocessing of endoscopes is estimated to be between 50-153 USD per cycle. (13)

The comprehensive reprocessing process is demonstrated in a simplified overview on page 9.

# Comprehensive reprocessing is complex, time consuming, and costly



# Recent cases and articles demonstrate a cross-contamination risk from 0.6 to 4.6%

## T.D. Waite et al. 2016 article. UK case (16)

*"Pseudo-outbreaks of Stenotrophomonas maltophilia on an intensive care unit in England"*

Two pseudo-outbreaks occurred due to contaminated reusable bronchoscopes, affecting a total of 18 patients.

### Conclusion

*"most notably the change to single-use bronchoscopes, have negated the false-positive reporting of S. maltophilia. In turn, this has reduced the risk of inappropriate antibiotic use and isolation of patients, and has increased patient safety."*

## M. Guy et al. 2016 article. France case (17)

*"Outbreak of pulmonary Pseudomonas aeruginosa and Stenotrophomonas maltophilia infections related to contaminated bronchoscope suction valves, Lyon, France, 2014"*

A total of 157 patients exposed to 216 bronchoscopic procedures from 1. December 2013 to 17. June 2014 were analysed. 10 cases of cross-contamination were linked directly to two bronchoscope suction valves, resulting in an **overall contamination risk of 4,6%**.

## J. Kovaleva et al. 2013 (2)

Review *"Transmission of Infection by Flexible Gastrointestinal endoscopy and Bronchoscopy"*

**Result:** Out of a total patient population of 569 the same contaminant was found in the patient as well as in the bronchoscope in 115 cases resulting in an **infection risk of approx. 20%**.

## C.J. Terjesen J. Kovaleva L. Ehlers 2017 (18)

*"Early Assessment of the Likely Cost Effectiveness of Single-Use Flexible Video Bronchoscopes"*

### Overall conclusions

- Using the current technology (**reusable bronchoscopes**) is estimated to have an average cost of **\$US424** and to hold a **0.7% risk of infection**. The newer technology (Single-Use) has an average cost per use of \$US305 and a 0% risk of infection.
- **Results show a possible saving of \$US118.56 per procedure and the elimination of a 0.7% risk of infection if the single-use option is adopted instead of the current technology.**

## Cori L. Ofstead et al. 2013 article. North America investigation (19)

*"Transmission of multidrug-resistant organisms and other pathogens via contaminated endoscopes: Can buildup biofilm be eliminated by routine cleaning and high-level disinfection?"*

**Result:** 251 bronchoscopes tested and in **4% of the bronchoscopes** organic debris remains after cleaning was found.

## Gavalda et al. 2015 article. Spain and Australia investigation (20)

*"Microbiological monitoring of flexible bronchoscopes after high-level disinfection and flushing channels with alcohol: Results and costs."*

**Result:** A total of **620 samples** were obtained. 56 samples (9%) tested positive for at least one specimen, of whom 3% were pathogenic or potentially pathogenic microorganisms. **Risk of contamination was 4,1%** without flushing channels with alcohol and **0,6%** when scope channels were flushed.

## Cori L. Ofstead et al. 2016 article. US case (21)

*"Practical toolkit for monitoring endoscope reprocessing effectiveness: Identification of viable bacteria on gastroscopes, colonoscopes, and bronchoscopes"*

**Result:** **4 out of 5 patient ready bronchoscopes tested positive for pathogen growth.**

# Sources

1. FDA MAUDE database <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfmaude/search.cfm>
2. J. Kovaleva et al. 2013; Review "Transmission of Infection by Flexible Gastrointestinal endoscopy and Bronchoscopy" Clinical Microbiology Revives April 2013 vol. 26 no. 2 231-254
3. ECRI institute <https://www.ecri.org>
4. US senate report January 13 2016 "Preventable Tragedies: Superbugs and How Ineffective Monitoring of Medical Device Safety Fails Patients" <https://www.help.senate.gov/imo/media/doc/Duodenoscope%20Investigation%20FINAL%20Report.pdf>
5. CDC Guideline 2008. Disinfection and Sterilization in Healthcare Facilities
6. FDA 2015 Materials of the Gastroenterology-Urology Devices Panel May 14-15, 2015: Meeting Materials Ph.D William Rutala presentation <https://www.fda.gov/advisorycommittees/committeesmeetingmaterials/medicaldevices/medicaldevicesadvisorycommittee/gastroenterology-urologydevicespanel/ucm445590.htm>
7. FDA Reprocessing of Reusable Medical Devices: Information for Manufacturers <https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/ReprocessingofReusableMedicalDevices/>
8. AAMI Releases 'Must-Have' Guide for Endoscope Reprocessing <http://www.aami.org/newsviews/newsdetail.aspx?ItemNumber=2468>
9. CDC safety alert September 11, 2015 "Immediate Need for Healthcare Facilities to Review Procedures for Cleaning, Disinfecting, and Sterilizing Reusable Medical Devices" <https://emergency.cdc.gov/han/hano0382.asp>
10. FDA safety communication September 17, 2015 "Infections Associated with Reprocessed Flexible Bronchoscopes: FDA Safety Communication" <https://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm462949.htm>
11. FDA safety communication November 13 2015 "FDA orders recall under consent decree for all Custom Ultrasonics automated endoscope reprocessors" <https://www.fda.gov/NewsEvents/Newsroom/PressAnnouncements/ucm472605.htm>
12. CDC release from the Healthcare Infection Control Practices Advisory Committee (HICPAC) January 25, 2017. "Essential Elements of a Reprocessing Program for Flexible Endoscopes – Recommendations of the HICPAC" <https://www.cdc.gov/hicpac/recommendations/flexible-endoscope-reprocessing.html>
13. Cori L. Ofstead et al 2017. "A GLIMPSE AT THE TRUE COST OF REPROCESSING ENDOSCOPES: RESULTS OF A PILOT PROJECT" In International Journal of Healthcare Central Service Material Management (www.iahcsmm.org)
14. Ofstead CL, Wetzler HP, Doyle EM, et al. "Persistent contamination on colonoscopes and gastroscopes detected by biologic cultures and rapid indicators despite reprocessing performed in accordance with guidelines. Am J Infect Control 2015, 43: 794–801.
15. FDA Guidance March 17 2015 "Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry and Food and Drug Administration Staff" <https://www.fda.gov/downloads/medicaldevices/deviceregulationandguidance/guidancedocuments/ucm253010.pdf>
16. T.D. Waite et al. 2016 article. UK case "Pseudo-outbreaks of Stenotrophomonas maltophilia on an intensive care unit in England" journal of hospital infection April 2016, Volume 92, Issue 4, Pages 392–396
17. M. Guy et al. "Outbreak of pulmonary Pseudomonas aeruginosa and Stenotrophomonas maltophilia infections related to contaminated bronchoscope suction valves" Europe's journal of infectious diseases epidemiology, prevention and control Eurosurveillance, Volume 21, Issue 28, 14 July 2016
18. C.J. Terhesen J. Kovaleva L. Ehlers "Early Assessment of the Likely Cost Effectiveness of Single-Use Flexible Video Bronchoscopes" PharmacoEconomics Open (2017) 1: 133. doi:10.1007/s41669-017-0012-9
19. Cori L. Ofstead al. "Transmission of multidrug-resistant organisms and other pathogens via contaminated endoscopes: Can buildup biofilm be eliminated by routine cleaning and high-level disinfection?" 2013, Ofstead & Associates
20. Gavalda et al. "Microbiological monitoring of flexible bronchoscopes after high-level disinfection and flushing channels with alcohol: Results and costs." Respiratory Medicine August 2015, Volume 109, Issue 8, Pages 1079–1085
21. Cori L. Ofstead et al. "Practical toolkit for monitoring endoscope reprocessing effectiveness: Identification of viable bacteria on gastroscopes, colonoscopes, and bronchoscopes" AJIC July 1, 2016, Volume 44, Issue 7, Pages 815–819
22. FDA, warning letters, <https://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/ReprocessingofReusableMedicalDevices/ucm454630.htm>
23. CDC safety Alert <https://emergency.cdc.gov/han/hano0382.asp>