



All-In-One - IHC, IF & SS
Clean Crisp Intense Stains
200 Slides in 8 hour Shift
Use your Preferred Antibodies



About i6000[™] ELITE DX

i6000™ Diagnostics System offers efficient cost effective solution with unsurpassed quality consistency, reliability, high capacity and throughput

Adaptable, Efficient, Productive and Profitable

- Dispense reagents as low as 100µl
- Adaptable to your work-flow-60 IHC slides in <3 hours
- Predictable cost structure
- High Capacity 60 slides and 53 antibodies per run
- More productivity per square foot of lab space
 - High throughput Up to 200 slides in a 8 hrs. shift
 - Start delayed run and return to stained slides next morning

Flexible System

- Multiple slide processing options Random, Continuous and STAT
- Multi-format specimen processing
 - FFPE or Frozen tissues
 - Cell preparations, Fine Needle Aspirates, Smears and more...
- Open system allows to use 3rd party antibodies

Easy to Use

- Advanced XYZ robotics with slide carriers and reagent racks allow quick set-up
- Top loading system ensures easy handling and management of slides and reagents
- Color coded GUI with real time assay parameter display for all slides
- Short learning curve allows cross-training and rotation of personnel
- Real time audio and visual alerts

Reproducible and Reliable Results

- Standardized staining process
- Over 300+ optimized protocols with ready to use reagents in barcoded vials
- Liquid level sensor ensures accurate reagent dispensing
- Disposable pipettes eliminate cross contamination

Reports and LIS connectivity

- Customized or standard reports for inventory management and regulatory compliance
- LIS connectivity*
 - · Capability to log sample remotely
 - Remote login of sample
 - Report access based on Sample ID or Request ID
- Step-by-step protocol log for every slide
- Run Parameters, reagent expiration dates and specimen used
- Reagent inventory management

The System



Removable Reagent Racks (2x30 vials)

• Offers quick and convenient storage of reagent vials Enables multiple staining protocols in a single run




Bar-coded Slides and Reagent Vials

- Requires minimal hands-on time
- Reduces data entry errors



Disposable Pipette Tips

• Eliminates cross-contamination of reagents




Removable Slide Racks with Splash Guard (5x12 Slides)

• Provides contamination free staining




Manage your resources and work flow efficiently

Special Stains




- A. Amyloid (Congo Red) Stain
- B. Colloidal Iron stain
- C. Alcian Blue/PAS Stain



- A. Mucicarmine Stain
- B. Gomori's Trichrome Blue Stain
- C. Reticulin Nuclear Fast Red Stain

Multiplex Stains




- A. SCC Anti-p63 + Anti-CK5
- B. Tonsil Anti-CD3 + Anti-CD20
- C. Metastatic Carcinoma in Lung Anti-CDX2 + Anti-CK7



Manage your resources and work flow efficiently

Multiplex Stains




A. anti-CD5 + anti-CD23


B. anti-CD4 + anti-CD8

C. anti-Ki67 + anti-Lambda

IHC



- A. Anti-CK18
- B. Anti-GCDFP-15
- C. Anti-ErbB-2/Her2



- A. MUC5AC in Gastro-intestine
- B. T-Cell in Tonsil
- C. Glypican-3 in Hepatocellular Ca

Integrated Workflow Solution







<u>**FZ**-Retriever</u>° **System**

Pre-treatment and Antigen Retrieval System



1) UU Plus[™] LIS Enabled Barcode Label Printer

- The i6000[™] is designed to work seamlessly with EZ-Retriever® providing Eco-friendly De-waxing, Rehydration and Antigen Retrieval
- Integrated Digital-data tracking system with i500 Plus™
- Broad range of optimized ready-to-use reagents

Ordering Information:

Name	Cat. No.
i6000 [™] ELITE Dx	AS6030
i6000 [™] ELITE Dx Software	4812-00093

Contact Information:

Regulatory Claims & Technical Support at: support@biogenex.com

USA:	49026 Milmont Drive, Fremont, CA 94538, USA		T: +1-800-421-4149 F: +1-510-824-1490
India:	SEZ Unit, Plot No.4, APIIC SEZ for Aerospace & Precision Engineering Adibatla, R. R. Dist 501510, Telangana State, India		T: +91-9676777714 F: +91-40-66529200
China:	Unit 511B, 5F Front Hall, Shanghai Exhibition Center No.1000 Middle Yan An Road, Shanghai 200041, PRC, China		T: +86-21-2220-6680 F: +86-21-2220-6679
Italy:	Via Savoldini, n.73, 24057, Martinengo, Bergamo, Italy	I	T: +39-036-398-6463



www.biogenex.com e-mail: customer.service@biogenex.com

X BioGenex_

© 2016 BioGenex Laboratories, Inc. All Rights Reserved