

## Logical®

### ORDERING INFORMATION

Product Code	Product Description	Units/ Case
MX9605	DPS SINGLE LINE PRESSURE MONITORING KIT 20/BX	20
MX9622	DPS DOUBLE LINE PRESSURE MONITORING KIT 20/BX	20
MX960P1	PRESSURE TRANSDUCER REUSABLE LOGICAL	1
MX960FR	LOGICAL DOME WITH TRIGGER FLUSH 25/BX	25
MX960XYP1	DPS DOME FEMALE/MALE 50/BX	50
MX261P1	CLAMP-DPS FOR BRACKET	1
MX262P1	BRACKET-DPS FOR 2 TRANSDUCERS	1
MX961Z50	LOGICAL@CABLE: WITHOUT CONNECTOR 1/EA	1

## Transtar®

### ORDERING INFORMATION

Product Code	Product Description	Units/ Case
MX9505T	TRANSTAR® 84" SINGLE LINE MONITORING KIT	10
MX9504T	TRANSTAR® 60" SINGLE LINE MONITORING KIT	10
MX9502T	TRANSTAR® DOUBLE LINE MONITORING KIT	10
MX950SCFT	TRANSTAR® DISPOSABLE PRESSURE TRANSDUCER WITH STOPCOCK AND 3CC/HR.FLUSH	10
MX800	MODULAR TRANSDUCER MOUNTING PLATE	20

## C-Fusor®

### ORDERING INFORMATION

Product Code	Product Description	Units/ Case
MX4805P1	C-FUSOR® 500 ML PRESSURE INFUSOR	1
MX4810P1	C-FUSOR® 1000 ML PRESSURE INFUSOR	1

## Clear-Cuff®

### ORDERING INFORMATION

Product Code	Product Description	Units/ Case
MX4705P1	CLEAR-CUFF® 500 ML PRESSURE INFUSOR	10
MX4710P1	CLEAR-CUFF® 1000 ML PRESSURE INFUSOR	10

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

#### Smiths Medical India Private Limited

508/509, Western Edge-2, W.E. Highway,  
Borivali (E), Mumbai - 66, India.  
T : + 91 22 28705210 F : +91 22 28705211  
E: csindia@smiths-medical.com  
www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group.

Smiths brand names Portex, and the Portex and Smiths Medical design marks are trademarks of Smiths Medical.  
Literature No: 0315/MDX/01/4P

smiths medical



smiths medical  
bringing technology to life

## Pressure Monitoring System

Transtar®



LogiCal®



Clear-Cuff®



C-Fusor®

For the use of a Registered Medical Practitioner or a Hospital or a Laboratory only.

medex®

## Pressure Transducer Systems

### LogiCal®

- LogiCal® transducer station remains in place as an extension to the monitor, while the disposable sets remain with the patient during transfers
- Built-in SURE-CAL® facility provides an integrity check of the system in seconds
- Disposable pressure cartridge with a smooth and narrow inner diameter ensures a bubble-free circuit
- Disposable patient line cartridge with integrated trigger flush and easy set up
- Guide rail with an audible click insures a correct dry coupling of the disposable pressure cartridge and the pressure transducer
- Unique design of the cable connector prevents fluid from entering the electronics parts
- Modular mounting system allows flexibility for multi-line monitoring kits



### Transtar®

- Ergonomic trigger flush device
- Transducer chip is integrated into the patient disposable
- Dependable and simple priming with straight fluid path
- Modular transducer mounting plate set-up
- Accurate to AAMI specifications and industry standards



## Pressure Infusor Devices

### C-Fusor®

- Crystal clear material allows maximum visibility of fluid levels
- Accurate and easy-to-read pressure gauge
- Allows for even pressure around the fluid container
- Three-way stopcock prevents air leaks
- Secure fastening
- Unique design ensures rapid set-up and fluid bag replacement
- Durable and easy to clean
- Large squeeze bulb for rapid and reliable pressure control
- 1000ml and 500ml sizes available



### Clear-Cuff®

- Clear - The thick and durable polyurethane construction provides good visibility of the fluid bag and fluid level
- Consistent Fluid Flow - For effective fluid management
- Cost-Effective - The Clear-Cuff® Pressure Infusor can be thoroughly cleaned and reused if desired (cleaning the Clear-Cuff® pressure infusor should be done according to the instructions for Use)
- Convenient - The built-in hook allows the fluid bag to be replaced without removing Clear-Cuff® from the I.V. pole
- Fast and Easy - The extra large squeeze bulb provides for fast inflation of the infusor
- The Hi-Flow™ venting stopcock allows Clear-Cuff® to be deflated almost instantly without squeezing the Clear-Cuff® bag

